


Product Description

Range : (00000.000 - 99999.999)Seconds.

EN-DTC has been specifically designed to calibrate Stopwatches , Timers and time stamp measurement.



Technical Specifications

Model No.	EN-DTC
Input Supply	230VAC/50Hz.
Parameter	Time
Counting Direction	UP
Range	(00000.000 - 99999.999)Seconds.
Resolution	0.001Second
Accuracy	± 0.01 % Of Reading.
Display	8 Digit LCD.
Out Put	Relay C/O With Rec. To Timing.

All product specifications are subject to change without prior notice.