

Product Description

Range : (Amb. ~ +250)°C



Applications

We offer best quality Low Range Dig. Wet Temperature Bath. Also known as EN-WTB250 , these are further more highly appreciated for portability and very very light weight. The rugged features of these products makes them highly sought after in the petrochemical , pharmaceutical industries & calibration laboratories etc. Our Baths are used for testing thermometers , sensors & thermocouples ranging from (Amb. ~ +250)°C. The total weight of these is 3.0 Kg.

Technical Specifications

Model No.	EN-WTB250
Temperature Range	(Amb. ~ +250)°C
Heating Specifications	
Liquid Tank Usable Inner Working Area	130mm. deep X 70mm. Well Dia.
Controller Specifications	
Type	Microprocessor Logic Controller
Display Range	(0 ~ 250)°C
Accuracy	±0.5°C
Resolution	0.1°C
Stability	±0.1°C
Uniformity	±0.2°C
Units	°C
Stabilization Time	Please allow 15 minutes to stabilize at particular Set Temp.
Displays	LED with SV & PV

All product specifications are subject to change without prior notice.

**Technical Specifications****Standard Specifications**

Heating Time	Time to reach 250°C from Ambient will take minimum 90 mins.
Switch Test Facility	Provided
Power Supply	230V AC \pm 5V
Dimensions	(270 X 140 X 200)mm.
Weight	3kgs.

All product specifications are subject to change without prior notice.